	MEM (Memo)	Dokukennung_Titel_Revision (Code_Titel_Revision) MTS_Treatment_Protocol_E_B	
	DG100 – Treatment protocol Revision 2013	Ausgabedatum (Date of Release) (TT.MM.JJ) 00.00.00	Page 1 / 8

Analysis of treatment protocols: application of ESWT for wound healing disorders

Shock waves are high-energy acoustic waves characterized by an initial rise of high pressure within a very short rise time, followed by a low tensile amplitude.

Extracorporeal shock wave therapy has been documented to have various effects on bones and soft tissue. ESWT acts through mechanotransduction and leads to biological responses at cellular level including the release of angiogenic growth factors known to play an important role in wound healing, associated with neovascularization, improved blood supply, cell proliferation and a significant bacterial effect.

Recent researches show that shock waves have highly beneficial effects on skin lesions since they stimulate the healing processes and in acute conditions appears to initiate a more rapid and effective healing phase.

Shockwave energy, frequency, number of pulses and number of re-treatments are important characteristics of treatment description in order to compare different ESWT studies and standardize shock wave treatment for various indications.

For shock waves application on wounds the number of pulses depends on the wound size.

Furthermore, the number of treatments is very case specific and depends on the type of wound healing disorders: some wounds react very well with fewer treatments while others require more treatments.

The following table shows an overview of recently published literature and MTS studies for wound indications, showing positive results and proving the effectiveness of ESWT in wound healing:

HUMAN STUDIES

MTS Spark Wave Therapy - Scientific Research: Wounds

RESEARCH INSTITUTION	INDICATION	PUBLICATION	AUTHOR	JOURNAL	IMPULSES AND ENERGY FLUX DENSITY	RESULTS
Centre for Severe Burns with Plastic Surgery, Unfallkrankenhaus Berlin, Berlin, Germany	Second degree burn wounds Clinical application	Prospective randomized phase II trial of accelerated reepithelialization of superficial second-degree burn wounds using extracorporeal shock wave therapy.	Ottomann C, Stojadinovic A, Lavin P, et al.	<i>Ann Surg.</i> 2012; 255(1):23-9.	100 impulses per cm² 0.1 mJ/mm²	<i>Pos.</i> 100%
Institute for Neuropalliative Rehabilitation, Royal Hospital for Neuro-disability, London, UK	Chronic decubitus ulceration Clinical application	Randomized control of extracorporeal shock wave therapy versus placebo for chronic decubitus ulceration.	Larking AM, Duport S, Clinton M, et al.	<i>Clin Rehabil.</i> 2010; 24(3):222-9.	200 + 100 impulses per cm² 5 pulses per sec 0.1 mJ/mm²	<i>Pos.</i> - Completely healed: 74% - 10% had 50% of epithelialization
Centre for Severe Burns with Plastic Surgery, Unfallkrankenhaus Berlin, Berlin, Germany	Re-epithelization of skin graft donor sites Clinical application	Prospective randomized trial of accelerated re-epithelization of skin graft donor sites using extracorporeal shock wave therapy.	Ottomann C, Hartmann B, Tyler j, et al.	<i>J Am Coll Surg</i> 2010; 211(3):361-7.	100 impulses per cm² 0.1 mJ/mm²	<i>Pos.</i>
Chang Gung Memorial Hospital, Kaohsiung Medical Center, Chang Gung University College of Medicine, Kaohsiung, Taiwan	Diabetic foot ulcers Clinical application	Extracorporeal shockwave treatment for chronic diabetic foot ulcers.	Wang C, Kuo Y, Wu R, et al.	<i>J Surg Res.</i> 2009; 152:96-103.	300 + 100 impulses per cm² 0.11 mJ/mm²	<i>Pos.</i> - Completely healed: 31% - Improved: 58%
Centre for Severe Burns with Plastic Surgery, Unfallkrankenhaus Berlin, Berlin,	Thermal lesion Clinical application	Accelerated reepitheliasation of a IIb° scald through extracorporeal shock wave therapy	Ottoman C, Thiele R, Hartmann B	<i>GMS Verbrennung medizin, 2009, 3: Doc01</i>	100 impulses per cm² 0.1-0.14 mJ/mm²	<i>Pos.</i>



MEM
(Memo)

Dokukennung_Titel_Revision
(Code_Titel_Revision)
MTS_Treatment_Protocol_E_B

**DG100 – Treatment protocol
Revision 2013**

Ausgabedatum
(Date of Release)
(TT.MM.JJ)
00.00.00

Page
3 / 8


Germany						
Department of Cardiothoracic Surgery, Medical University of Vienna, Austria	Wound healing after vein harvesting Clinical application	Prophylactic low-energy shock wave therapy improves wound healing after vein harvesting for coronary artery bypass graft surgery: A prospective, randomized trial.	Dumfarth J, Zimpfer D, Vögele-Kadletz, et al.	<i>Ann Thorac Surg.</i> 2008; 86:1909-13.	25 per cm 0.1 mJ/mm²	<i>Pos.</i>
AUVA-Trauma Center Meidling, Vienna, Austria	Chronic wounds Clinical application	Shock wave therapy for acute and chronic soft tissue wounds: A feasibility study.	Schaden W, Thiele R, Kolpl C, et al.	<i>J Surg Res.</i> 2007; 143:1-12.	100 per cm² 0.1 mJ/mm²	<i>Pos.</i> 75%

ANIMAL STUDIES

MTS Spark Wave Therapy - Scientific Research: Wounds


RESEARCH INSTITUTION	INDICATION	PUBLICATION	AUTHOR	JOURNAL	IMPULSES AND ENERGY FLUX DENSITY	RESULTS
Ludwig Boltzmann Institute for Experimental and Clinical Traumatology, Vienna, Austria	Tissue revascularisation Animal model	Extracorporeal shockwave therapy (ESWT) minimizes ischemic tissue necrosis irrespective of application time and promotes tissue revascularization by stimulating angiogenesis.	Mittermayr R, Hartinger J, Antonic V, et al.	<i>Ann Surg.</i> 2011; 253(5):1024-32.	300 impulses per cm² 5 pulses per sec 0.1 mJ/mm²	<i>Pos.</i>
Chang Gung Memorial Hospital, Kaohsiung Medical Center, Chang Gung University College of Medicine,	Enhanced wound healing in rat model of STZ-induced diabetes Animal model	Extracorporeal shock-wave therapy enhanced wound healing via increasing topical blood perfusion and tissue regeneration in a rat model of	Kuo, YR, Wang CT, Wang FS et al.	<i>Wound Rep Reg.</i> 2009; 17:522-530.	800 impulses per cm² 0.9 mJ/mm²	<i>Pos.</i>

Kaohsiung, Taiwan		STZ-induced diabetes.				
Naval Medical Research Center, Combat Casualty Care Department of Regenerative Medicine, Silver Spring, USA	Burn injury Animal model	Extracorporeal shock wave therapy suppresses the early proinflammatory immune response to a severe cutaneous burn injury.	Davis T, Stojadinovic A, Anam K, et al.	<i>Int Wound J.</i> 2009; 6:11-21.	200 impulses per cm² 5 pulses per sec 0-1 mJ/mm²	<i>Pos.</i> - Significant reduction of neutrophils : 60-68% - Significant reduction of macrophages: 55-66%
Walter Reed Army Medical Center, Combat Wound Initiative Program, Department of Surgery, Washington, USA	Angiogenic response in skin isografts Animal model	Angiogenic response to extracorporeal shock wave treatment in murine skin isografts.	Stojadinovic A, Elster E, Anam EA, et al.	<i>Angiogenesis.</i> 2008; 11(4):369-380.	200 impulses per cm² 3 pulses per sec 0.1 mJ/mm²	<i>Pos.</i>
University of Kentucky Chandler Medical Center, Division of plastic and reconstructive surgery, Department of Internal Medicine and College of Dentistry, Lexington, KY	Skin flap model Animal model	A study of the biological factors and wound healing of a skin flap model treated with unfocused extracorporeal shockwave therapy	Edelmann S, Fink BF, et al.	<i>University of Kentucky Chandler Medical Center</i>	500 impulses per cm² 0.15 mJ/mm²	<i>Pos.</i>
Ludwig Boltzmann Institute for experimental and clinical Traumatology, Austrian Cluster for Tissue Regeneration –Trauma Center Meidling, AUVA, Vienna, Austria	Skin flap model Animal model	How many shockwaves are enough? Dose-response relationship in ischemic challenged tissue	Mittermayr R, Hartinger J, Hofmann M, et al.	<i>Annual Congress ISMST 2008, Abstract No. 64</i>	30 impulses per animal 300 impulses per animal 1000 impulses per animal	<i>Best treatment:</i> 300 impulses per animal

	MEM (Memo)	Dokukennung_Titel_Revision (Code_Titel_Revision) MTS_Treatment_Protocol_E_B	
	DG100 – Treatment protocol Revision 2013	Ausgabedatum (Date of Release) (TT.MM.JJ) 00.00.00	Page 5 / 8

Comparison treatment protocols Wang et al. (energy-impulses-intervals)










	Wang et al. (2009) “Extracorporeal shock wave treatment for chronic diabetic foot ulcers”	Wang et al. (2011) “Treatment of diabetic foot ulcers: a comparative study of extracorporeal shockwave therapy and hyperbaric oxygen therapy”
	MTS Orthowave180 New: Orthogold100	Sanuwave DermaPACE
Indication	Chronic diabetic foot ulcers	Chronic diabetic foot ulcers
Impulses	300 + (cm² x 10) Total surface	500 + (cm² x 8) 1 cm plus in all directions
Energy Flux Density	Unfocused 0.11 mJ/mm²	Focused 0.23 mJ/mm²
Intervals	Once every 2 weeks (total 3 treatments in 6 weeks)	2 times per week for 3 weeks (total of 6 treatments)


	MEM (Memo)	Dokumennummer_Titel_Revision (Code_Titel_Revision) MTS_Treatment_Protocol_E_B	
	DG100 – Treatment protocol Revision 2013	Ausgabedatum (Date of Release) (TT.MM.JJ) 00.00.00	Page 6 / 8

PROTOCOL CALCULATION ANALYSIS

		MTS Dermagold100	Sanuwave DermaPACE	X wound size: (1x2) cm ²	Y wound size: (2x5) cm ²	Z wound size: (4x5) cm ²
O L D	MTS A) IFU	100 shocks/cm ²		157	785	1587
	MTS B) User guide	350 + (cm ² x 10)		365,7	428	507
N E W	MTS C) Proposal I	50 shocks/cm ² Total surface		100	500	1000
	MTS D) Proposal II	350 + (cm ² x 30) Total surface		410	650	950
	MTS E) Proposal III	30 shocks/cm ² Total surface		60	300	600
	MTS F) Proposal IV	350 + (cm ² x 20) Total surface		390	550	750
	Sanuwave G) Measurement		500 + (cm ² x 8) 1 cm plus in all directions	596	724	836

Examples of calculation impulses

	MTS old	MTS new	Sanuwave														
X	 Wound size (1x2) cm ² total surface 1.57 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">A</td><td>157</td></tr> <tr><td style="background-color: #cccccc;">B</td><td>365,7</td></tr> </table>	A	157	B	365,7	 Wound size (1x2) cm ² total surface 2 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">C</td><td>100</td><td style="background-color: #cccccc;">E</td><td>60</td></tr> <tr><td style="background-color: #cccccc;">D</td><td>410</td><td style="background-color: #cccccc;">F</td><td>390</td></tr> </table>	C	100	E	60	D	410	F	390	 Wound size (1x2) + 1 cm in all directions: 12 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">G</td><td>596</td></tr> </table>	G	596
A	157																
B	365,7																
C	100	E	60														
D	410	F	390														
G	596																
Y	 Wound size (2x5) cm ² total surface 7.85 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">A</td><td>785</td></tr> <tr><td style="background-color: #cccccc;">B</td><td>428</td></tr> </table>	A	785	B	428	 Wound size (2x5) cm ² total surface 10 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">C</td><td>500</td><td style="background-color: #cccccc;">E</td><td>300</td></tr> <tr><td style="background-color: #cccccc;">D</td><td>650</td><td style="background-color: #cccccc;">F</td><td>550</td></tr> </table>	C	500	E	300	D	650	F	550	 Wound size (2x5) + 1 cm in all directions: 28 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">G</td><td>724</td></tr> </table>	G	724
A	785																
B	428																
C	500	E	300														
D	650	F	550														
G	724																
Z	 Wound size (4x5) cm ² total surface 15.7 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">A</td><td>1587</td></tr> <tr><td style="background-color: #cccccc;">B</td><td>507</td></tr> </table>	A	1587	B	507	 Wound size (4x5) cm ² total surface 20 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">C</td><td>1000</td><td style="background-color: #cccccc;">E</td><td>600</td></tr> <tr><td style="background-color: #cccccc;">D</td><td>950</td><td style="background-color: #cccccc;">F</td><td>750</td></tr> </table>	C	1000	E	600	D	950	F	750	 Wound size (4x5) + 1 cm in all directions: 42 cm ² <table border="1" style="margin-top: 5px;"> <tr><td style="background-color: #cccccc;">G</td><td>836</td></tr> </table>	G	836
A	1587																
B	507																
C	1000	E	600														
D	950	F	750														
G	836																

	MEM (Memo)	Dokukennung_Titel_Revision (Code_Titel_Revision) MTS_Treatment_Protocol_E_B	
	DG100 – Treatment protocol Revision 2013	Ausgabedatum (Date of Release) (TT.MM.JJ) 00.00.00	Page 8 / 8

Conclusions:

The different way to calculate the surface of the wound and the number of required impulses can lead to different results!

As reported in the analysed literature studies (both human and animal studies), most of them showed a number of impulses equal to 100 impulses per cm².

There is only one scientific evidence for a doses-outcome-relation from Mittermayr et al., who investigated skin flap outcome in respond to various total amounts of impulses, in an animal flap model.

In the ischemic area of a rodent epigastric flap, different amounts of total shock wave impulses were applied (30, 300, and 1000) to reduce tissue necrosis. They found out that the best dose is 300 impulses. Thus, they showed that more impulses are not necessary and that a higher amount does not improve the final outcome.

The proposed protocol reflects also the latest Wang et al. study (2011) in which it is shown that the protocol with 8 impulses per cm² is similar to the one used for Dermapace, but the energy flux density is different: in Wang study they used the double amount of energy.

Based on this scientific paper, a certain minimum dose is required and the maximum doses should be limited. Therefore it is justified to limit the impulses with 20 impulses per cm², plus a minimum quantity equal to 350 shocks (Proposal IV). This means to have more shocks on small wounds, and less shocks in large wounds compared to a linear formula, used in previous protocols.

Konstanz, 20.06.2013
Nicoletta Fotino / Ralph Reitmajer